

Effects of Direct Oral Anticoagulants on Hemostasis testing

About this guide

The ISTH SCC¹ and British Committee for Standards in Haematology² recommend that laboratories understand the effects of direct oral anticoagulants (DOACs) on Hemostasis assays. This guide provides detailed reference data to illustrate the impact of DOACs on HemosIL[®] assays.

What has changed in anticoagulant therapy?

Hemostasis testing is critical in the overall clinical assessment of patients requiring anticoagulant therapy. Warfarin has been the gold standard in oral anticoagulant therapy for many years. DOACs, such as dabigatran, rivaroxaban, apixaban, edoxaban and betrixaban,* have been approved for clinical use in many countries and do not require routine monitoring. Therefore, understanding the characteristics of DOACs and differentiating their effects on routine Hemostasis tests is essential.

DOACs vs. Warfarin^{3,8}

	Dosing	Time to Peak Anticoagulation	Half-life	Antidote	Renal Excretion (%)	Drug Interaction Potential
Warfarin	Once Daily	3-5 Days	20-60 Hours	Vitamin K	0	Multiple, Especially CP2C9 Inhibitors
Dabigatran Etexilate	Twice Daily	~2 Hours	12-17 Hours	Idarucizumab	~80	P-Glycoprotein Inhibitors and Inducers
Rivaroxaban	Once Daily	2-4 Hours	5-9 Hours	Andexanet alfa	~36	CYP3A4 and P-Glycoprotein, Inhibitors and Inducers
Apixaban	Twice Daily	3-4 Hours	12 Hours	Andexanet alfa	~25	CYP3A4 and P-Glycoprotein, Inhibitors and Inducers
Edoxaban	Once Daily	1-2 Hours	10-14 Hours	In Development	~50	CYP3A4 and P-Glycoprotein, Inhibitors and Inducers
Betrixaban	Once Daily	3-4 Hours	19-27 Hours	In Development	~11	P-Glycoprotein Inhibitors

* Betrixaban was recently FDA 510(k)-cleared and is not be available in all territories. There is currently limited data available on the effect of betrixaban on Hemostasis assays.

HOW do DOACs affect Routine Hemostasis tests?

PT, APTT, Fibrinogen and Thrombin Time (TT) results can be affected by DOACs. Potential effects are described below.

Effects of DOACs on PT, APTT, Fibrinogen, TT and D-Dimer Results

	PT	APTT	Fibrinogen Activity	TT	D-Dimer
Dabigatran Etexilate	Prolonged + ^{3†}	Prolonged ++ ^{3†}	No Effect or Falsely Low ^{6,#}	Prolonged +++ ³	No Effect ⁴
Rivaroxaban	No Effect or Prolonged ++ ^{1,3†}	No Effect or Prolonged ++ ^{1,3†}	No Effect ³	No Effect ³	No Effect ⁴
Apixaban	No Effect or Prolonged ++ ^{1,3†}	No Effect or Prolonged ++ ^{1,3†}	No Effect ³	No Effect ³	No Effect ⁴
Edoxaban	No Effect or Prolonged ++ ^{11†}	No Effect or Prolonged ++ ^{11†}	No Effect ¹¹	No Effect ¹¹	No Effect ¹¹

+ slight increase ++ moderate increase +++ marked increase

** Data on file, IL.

† Direct IIa/Xa inhibitors may show variable effects, depending on the drug and drug concentration. In addition, different reagents may have different sensitivities.

Effect is assay-method-and drug-dependent. Most fibrinogen assays will show no effect with dabigatran.

References

1. Baglin T, *et al.* Measuring oral direct inhibitors (ODIs) of thrombin and factor Xa: a recommendation from the Subcommittee on Control of Anticoagulation of the Scientific and Standardisation Committee of the International Society on Thrombosis and Haemostasis. *J Thromb Haemost.* 2013;11(4):756-60.
2. Baglin T, *et al.* Effects on routine coagulation screens and assessment of anticoagulant intensity in patients taking oral dabigatran or rivaroxaban: guidance from the British Committee for Standards in Haematology. *Br J Haematol.* 2012;159(4):427-9.
3. Funk DM. Coagulation assays and anticoagulant monitoring. *Hematology Am Soc Hematol Educ Program.* 2012;2012:460-5.
4. Favaloro EJ, *et al.* Interference of direct oral anticoagulants in haemostasis assays; high potential for diagnostic false positive and false negatives. *Blood Transfus.* 2017;7:1-4.
5. Mani H. Interpretation of coagulation test results under direct oral anticoagulants. *Int Jnl Lab Hem.* 2014;36:261-8.
6. Dager WE, *et al.* Dabigatran effects on the international normalized ratio, activated partial thromboplastin time, thrombin time, and fibrinogen: a multicenter, in vitro study. *Ann Pharmacother.* 2012;46(12):1627-36.
7. Hoxa A, *et al.* Detection of lupus anticoagulant in the era of direct oral anticoagulants. *Autoimmun Rev.* 2017;16(2):173-8.
8. Kim YA, *et al.* The effects of dabigatran on lupus anticoagulant, diluted plasma thrombin time, and other specialized coagulation assays. *Int J Lab Hematol.* 2015;37(4):e81-4.
9. Maryamchik E, *et al.* Rivaroxaban causes missed diagnosis of protein S deficiency but not of activated protein C resistance (factor V Leiden). *Arch Pathol Lab Med.* 2018; 142(1):70-4.
10. Gosselin R, *et al.* Comparison of the effect of the anti-Xa direct oral anticoagulants apixaban, edoxaban, and rivaroxaban on coagulation assays. *Int Jnl Lab Hem.* 2016;38:505-13.
11. Douxfils J, *et al.* Edoxaban: Impact on routine and specific coagulation assays. A practical laboratory guide. *Thromb Haemost.* 2015;115(2):368-81.

HEMOSTASIS INNOVATION IS HERE ►

The Instrumentation Laboratory logo, ACL, ACL AcuStar, ACL Elite, ACL TOP, HemoCell, HemoHub, HemosIL, ReadiPlasTin, RecombiPlasTin, SynthAFax and SynthASil are trademarks of Instrumentation Laboratory Company and/or one of its subsidiaries or parent companies and may be registered in the United States Patent and Trademark Office and in other jurisdictions. All other product names, company names, marks, logos and symbols are trademarks of their respective owners.
©2019 Instrumentation Laboratory. All rights reserved.

For more information, contact your local Werfen sales representative or distributor.

How do DOACs affect HemosIL PT and APTT assays?

Effect of DOACs on HemosIL PT Results

Drug Concentration (ng/mL)	Dabigatran			Rivaroxaban			Apixaban			Edoxaban		
	RdPTn*	RPT 2G	PT HS Plus	RdPTn*	RPT 2G	PT HS Plus	RdPTn*	RPT 2G	PT HS Plus	RdPTn*	RPT 2G	PT HS Plus
0	1.00	1.00	1.00	1.00	1.00	1.00	1.00	1.00	1.00	1.00	1.00	1.00
50	1.04	1.02	1.05	1.16	1.13	1.19	1.07	1.06	1.07	1.04	1.04	1.14
100	1.05	1.02	1.08	1.30	1.28	1.33	1.15	1.13	1.12	1.10	1.10	1.34
200	1.08	1.02	1.13	1.62	1.58	1.70	1.28	1.27	1.23	1.28	1.27	1.64
300	1.11	1.04	1.18	1.94	1.89	1.98	1.45	1.40	1.32	1.43	1.41	1.97
400	1.26	1.18	1.34	2.24	2.18	2.31	1.61	1.56	1.44	1.62	1.59	2.24
500	1.63	1.60	1.68	2.59	2.50	2.56	1.78	1.71	1.54	1.76	1.73	2.50

* Currently available on ACL TOP® Family and ACL TOP Family 50 Series Hemostasis Testing Systems.
 RdPTn = ReadiPlasTin®; RPT 2G = RecombiPlasTin 2G®
 PT ratios (drug concentration vs. no drug) measure the effect of the DOACs at increasing concentrations. 1.00 = no effect. Normal range PT ratio < 1.25.


Effect of DOACs on HemosIL APTT Results

Drug Concentration (ng/mL)	Dabigatran			Rivaroxaban			Apixaban			Edoxaban		
	SynthASil®	APTT-SP	SynthAFax®	SynthASil	APTT-SP	SynthAFax	SynthASil	APTT-SP	SynthAFax	SynthASil	APTT-SP	SynthAFax
0	1.00	1.00	1.00	1.00	1.00	1.00	1.00	1.00	1.00	1.00	1.00	1.00
50	1.41	1.47	1.39	1.18	1.14	1.17	1.09	1.09	1.08	1.17	1.07	1.13
100	1.64	1.66	1.70	1.26	1.22	1.29	1.14	1.14	1.13	1.22	1.15	1.23
200	2.09	2.01	2.14	1.39	1.34	1.45	1.21	1.20	1.20	1.32	1.27	1.40
300	2.42	2.23	2.52	1.52	1.51	1.60	1.25	1.23	1.25	1.41	1.39	1.54
400	2.65	2.46	2.80	1.63	1.59	1.76	1.29	1.28	1.29	1.49	1.50	1.70
500	2.93	2.66	3.07	1.72	1.71	1.89	1.32	1.31	1.34	1.56	1.59	1.83

APTT ratios (drug concentration vs. no drug) measure the effect of the DOACs at increasing concentrations. 1.00 = no effect. Normal range APTT ratio < 1.25.

Highlighted areas show ratios above the normal range.


DOAC Effects on HemosIL ReadiPlasTin*



Legend: Dabigatran (blue square), Rivaroxaban (yellow triangle), Apixaban (purple circle), Edoxaban (teal square), Normal Range (orange line)

*Available on ACL TOP Family and ACL TOP Family 50 Series Hemostasis Testing Systems.

DOAC Effects on HemosIL SynthASil



Legend: Dabigatran (blue square), Rivaroxaban (yellow triangle), Apixaban (purple circle), Edoxaban (teal square), Normal Range (orange line)

DOACs target Factor IIa and Factor Xa, key proteins in the coagulation cascade that can interfere with Hemostasis specialty tests.

Effects of DOACs on Specialty Tests

	Dabigatran	Rivaroxaban	Apixaban	Edoxaban
Protein S (Antigen)	No Effect ⁵	No Effect ³	No Effect ³	No Effect ¹¹
Protein S (Clotting)	Falsely Increased ³	Falsely Increased ³	Falsely Increased ³	Falsely Increased ¹¹
Protein C (Chromogenic)	No Effect ³	No Effect ³	No Effect ³	No Effect ¹¹
Protein C (Clotting)	Falsely Increased ³	Falsely Increased ³	Falsely Increased ³	Falsely Increased ¹¹
Antithrombin (Xa)	No Effect ³	Falsely Increased ³	Falsely Increased ³	Falsely Increased ¹¹
Antithrombin (IIa)	Falsely Increased ³	No Effect ³	No Effect ³	No Effect ¹¹
Liquid Anti-Xa	No Effect ⁴	Increased ⁴	Increased ⁴	Increased ¹¹
dRVVT Screen & Confirm	Likely to Misclassify as LA ^{8†}	Possible to Misclassify as LA ^{7†}	Possible to Misclassify as LA ^{7†}	Possible to Misclassify as LA ^{11†}
PT-Based Factor Assays (FII, FV, FVII, FX)	Falsely Low ^{3,11}	Falsely Low ^{3,11}	Falsely Low ^{3,11}	Falsely Low ^{11,12}
APTT-Based Factor Assays (FVIII, FIX, FXI, FXII)	Falsely Low ^{3,4,11}	Falsely Low ^{3,4,11}	Falsely Low ^{3,4,11}	Falsely Low ^{11,12}
APC Resistance	Falsely Elevated Ratio ³	Mild False Elevations in Ratio, but Assay Remains Accurate ⁹	Mild False Elevations in Ratio, but Assay Remains Accurate ³	Mild False Elevations in Ratio, but Assay Remains Accurate ¹¹
Factor VIII (Chromogenic)	No Effect ³	Falsely Low ³	Falsely Low ³	Falsely Low ³
Factor XIII (Immunologic)	No Effect ⁵	No Effect ⁵	No Effect ¹⁰	No Effect ¹⁰

† Data on file, IL.

¹¹ Falsely low, may yield false impression of an inhibitor. LA = Lupus anticoagulant